

Examination Table Papers

Item Codes

Item Code	Description
80201	Smooth White
80202	Crepe White
80203	Smooth White
80204	Crepe White



Intended Use

Medical examination table paper provides protection with a one-sided glaze finish to repel moisture. All items are disposable and meant for single-use with each application.

Style

See details in table above

Material

Item 80201 & 80203 consist of pre-consumer secondary fiber. Not made with natural latex rubber
 Item 80202 & 80204 consist of Virgin Kraft hardwood/softwood blend and pre-consumer secondary fiber. Not made with natural latex rubber

Measuring

Item	Style	Width	Length	Center Core Diameter
80201	Smooth	21" (53 cm)	225 ft (68.5 m)	1 5/8" (4 cm)
80202	Crepe	21" (53 cm)	125 ft (38 m)	1 5/8" (4 cm)
80203	Smooth	18" (46 cm)	225 ft (68.5 m)	1 5/8" (4 cm)
80204	Crepe	18" (46 cm)	125 ft (38 m)	1 5/8" (4 cm)

Packaging

To withstand transportation and storage. Case made of suitable corrugate fiberboard.

Basis Weight

Item 80201 & 80203: (25 g/m²)
 Item 80202 & 80204: (30 g/m²)

Count

12 rolls per case

Country of Origin

USA

Storage Conditions

Store in a cool dry place