

ELASTIC BANDAGES

PRODUCT CATALOG

Our elastic bandages are designed to treat muscle strains and sprains by compression of the injured area. The bandages are available with either secure clip fasteners or hook-and-loop (Velcro) closures that provide an even and secure adherence to itself.

The elastic bandages are made from polyester without any natural rubber latex, come in both non-sterile and sterile variants, and allow for easy application and removal while maintaining elasticity for reapplication. Ideal for orthopedic and athletic use, as well as in the treatment of edemas and similar conditions, these bandages are an essential component of any medical inventory.



Training on the use of these bandages can be provided upon request and is available via many formats including in-person sessions, video media, web-based material, and written instructions.



Elastic Bandages with Clip and Hook-and-Loop (Velcro) Fasteners

ELASTIC BANDAGES

PRODUCT CATALOG

Product Portfolio and Specifications

Full Item Description	Part Number	Packaging
BANDAGE ELASTIC WITH CLIP FASTENER NON-STERILE 6"X 5YDS	B053642840ON	12 EA/PK, 20 PK/CS
BANDAGE ELASTIC WITH SELF-GRIPPING HOOK-AND-LOOP FASTENER NON-STERILE 6"X 5YDS	B053652840ON	12 EA/PK, 20 PK/CS
BANDAGE ELASTIC WITH CLIP FASTENER NON-STERILE 4"X 5YDS	B053642820ON	12 EA/PK, 30 PK/CS
BANDAGE ELASTIC WITH SELF-GRIPPING HOOK-AND-LOOP FASTENER NON-STERILE 4"X 5YDS	B053652820ON	12 EA/PK, 30 PK/CS
BANDAGE ELASTIC WITH CLIP FASTENER NON-STERILE 3"X 5YDS	B053642810ON	12 EA/PK, 40 PK/CS
BANDAGE ELASTIC WITH SELF-GRIPPING HOOK-AND-LOOP FASTENER NON-STERILE 3"X 5YDS	B053652810ON	12 EA/PK, 40 PK/CS
BANDAGE ELASTIC WITH CLIP FASTENER NON-STERILE 2"X 5YDS	B053642800ON	12 EA/PK, 60 PK/CS
BANDAGE ELASTIC WITH SELF-GRIPPING HOOK-AND-LOOP FASTENER NON-STERILE 2"X 5YDS	B053652800ON	12 EA/PK, 60 PK/CS
BANDAGE ELASTIC WITH SELF-GRIPPING HOOK-AND-LOOP FASTENER STERILE 4"X 5YDS	B533652820PH	1 EA/PK, 200 PK/CS
BANDAGE ELASTIC WITH SELF-GRIPPING HOOK-AND-LOOP FASTENER STERILE 6"X 5YDS	B533652840PH	1 EA/PK, 150 PK/CS
BANDAGE ELASTIC WITH CLIP FASTENER STERILE 4"X 5YDS	B533642820PH	1 EA/PK, 200 PK/CS
BANDAGE ELASTIC WITH CLIP FASTENER STERILE 6"X 5YDS	B533642840PH	1 EA/PK, 150 PK/CS
BANDAGE ELASTIC WITH SELF-GRIPPING HOOK-AND-LOOP FASTENER STERILE 2"X 5YDS	B533652800PH	1 EA/PK, 400 PK/CS
BANDAGE ELASTIC WITH CLIP FASTENER STERILE 3"X 5YDS	B533642810PH	1 EA/PK, 300 PK/CS
BANDAGE ELASTIC WITH CLIP FASTENER STERILE 2"X 5YDS	B533642800PH	1 EA/PK, 400 PK/CS